

LCA 30

The Leyro LCA 30 multifunction temperature calibrator is an equipment suitable for all conditions and applications, covering different kind of industries. The Leyro LCA 30 multifunction is the most advanced equipment in its class. Designed to work in a few minutes, as a dry block and as a liquid bath. Thanks to the wide variety of inserts available, it allows you to calibrate any probe or thermocouple. The LCA 30 device has a temperature range of -35°C ... 165°C *, stability of 0.01°C and accuracy of 0.1°C.



APPLICATIONS

Calibration Laboratories
Testing Labs
Universities / R & D
Chemical / Pharmaceutical Industry
Food industry
Maintenance and Service

HIGHLIGHTS

Fast Response Time
Temperature range -35 ... 165°C*
External Reference Temperature
Inmersion Depth: 190 mm
Resolution: 0,1/0,01/ 0,001°C
Stability: 0,01°C
Accuracy: 0,1°C
Uniformity: 0,05°C
Portable
3 Years Warranty

General Information

User interface with a touch screen of 4,3 inches All in one equipment (Dry Block/Liquid Bath) Data Register in a 80 Mb internal memory Function of 40 programmes

Trend Graph

Ethernet communication and two USB ports

Barcode reader

External Reference PT 100 4 Wires probe





Technical Data

Operating Temperature Range(°C) -35...165°C*
Accuracy 0,1°C

Resolution 0.1/0,01/0,001°C
Temperature Control Internal/ External
PT100 (3/4 Wires)

Calibration Volume Ø 60mm x 190 mm

Warm-Up time From -5 to 100° C | 25 Minutes Cooling time From 25 to -25° C | 45 Minutes

Housing and assembly Portable

Dimensions W x L x H (cm) 280 x 370 x 440

Weight (kg) 15.2 kg

Supply 100 ... 234 Vac 50/60 Hz

Power Consumption Approx 310 W

LIQUID BATH** DRY BLOCK

Temperature Stability ± 0.01 °C Temperature Stability 0,01 °C Uniformity ± 0.05 °C / 0,09 °F Uniformity 0,05 °C

Accuracy 0,1°C / 02°F Insert Dimensions 59,5mm x165mm

Well capacity 0.7 liters Standar Insert A. 1 x Ø 3,5 /1 x Ø 6,5/ 1 x Ø 8,5/ 1 x Ø 10,5

Optional Inserts

B. 2 x Ø 3,5 /2x Ø 4,5/ 2 x Ø 6,5/ 2 x Ø 8,5/ 2 x Ø 10,5

C. 3 x Ø 3,5 /3x Ø 6,5/ 3 x Ø 8,5/3 x Ø 10,5

 $D.\ 2\ x\ \varnothing\ 3,5\ /1\ x\ \varnothing\ 4,5/\ 1\ x\ \varnothing\ 5,0/1\ x\ \varnothing\ 5,5/1\ x\ \varnothing\ 6,5/\ 1\ x\ \varnothing\ 8,5\ /1\ x\ \varnothing\ 9,0/\ 1\ x\ \varnothing\ 9,5/1\ x\ \varnothing\ 10,5$


E. Blank

*Minimum Temperature Range is 55 °C, below ambient. Absolute minimum -35°C

**It is recommended to use 10 Cst silicone oil







Options





Scope of delivery

Basic Instrument LCA 30 Dry Block / Liquid Bath Function

Power Cord Trolley transport case with wheels and handle

Teflon Magnetic Stirrer ENAC Calibration Certificate

Manufacturing Traceable Calibration Certificate 8 Points Methacrylate lid with 6 inserts

Instruction Manual External Control
Drain Pump Pt 100/ 4 Wires

3 Years of Warranty against factory defects Several types of inserts (See the Dry Block Tecnical Data)